

Biohazardous Waste

Solids

Non-sharps

- dishes
- flasks
- tubes
- serological pipettes
- plastics
- pipette tips
- paper towels
- gloves



Discard into biohazard bag-lined garbage bin (IN BSC if BSL2)

Sharps

- broken glass
- glass vials
- metal caps
- needles
- syringes
- razors/scalpels
- Pasteur pipettes



Discard into puncture-proof sharps container (IN BSC if BSL2)

Discard **contained** solids in double lined biohazard boxes in waste room.

Liquids

Aspirate all liquids from solids before disposal.



Decontaminate liquids before emptying into the sink.

*** 30 minute minimum contact time with 10% bleach followed by 3 minutes running tap water



Plastic-ware is emptied of liquids and discarded into the biohazard bag-LINED BSC garbage container. (and gloves and paper towels)



Sharps, glass and metal are emptied of liquids and discarded into the BSC sharps container.



CONTAINMENT OF BIOHAZARDOUS SOLIDS





TIED in BSC



Spray with 70% ethanol




$\frac{3}{4}$ full MAX!!!



SECURELY SEALED in BSC


$\frac{3}{4}$ is FULL. Close inner bag.
Close outer bag and box
flaps. Return to waste room.

DISPOSAL OF **CONTAINED** **BIOHAZARDOUS SOLIDS**



BLEACH to 10% (for minimum of 30 minutes) or other biohazard - neutralizing chemical as per manufacturer's instructions

MIX WELL



Run tap 3 minutes !



DISPOSAL OF BIOHAZARDOUS LIQUIDS